PRODUCT DESCRIPTION

GEFRIERGETROCKNETE PRODUKTE ■ / FREEZE-DRIED PRODUCTS ■

1. Kortikale Streifen

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
Cortical Struts

Example 1: ALO 103 Cortical Struts 60 mm x 20 mm x > 2 mm

Corticalis struts produced from femur diaphysis, is the most callous part of the bone. Struts are used where immediate stability is required. Examples of use are the splinting of periprothetic fractures or augmentations of stress shielding zones.

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
Cortical Pins

Example 2: ALO 204 Cortical Pin 20 mm x 5 mm x 5 mm

Pins made of corticalis are applied as fixation of osteosynthesis as by discretion of the surgeon or if metal implants are not suited or absorbable implants provide too few stability. Advantage: slow incorporation (about 4 months) during the consolidation of osteosynthesis, no necessity of removal.

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


Example 4: ALO 311 Corticocancellous Granula 30 cc, 2 - 5 mm

Example 5: ALO 314 Cancellous Cubes 20 cc, 5 x 5 x 5 mm

Cancellous Granula is perfect for the filling of irregular cavitory defects. Examples therefore are being cysts as well as defects after loosening of endoprosthesis. The particle sizes are chosen according to the purpose of application, the spongiosa structure is always maintained and grants a quick incorporation. The application of **Corticocancellous Granula** is the same as for Cancellous Granula, but due to the 1:1 ratio cancellous cortical granula increased the stability is provided after impaction.

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
Primary anchorage at good incorporation.

Structural grafts are - alone or in combination with granula - suited for achieving partial load bearing stability.

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MEDICAL REFERENCES

Special preparations or enquiries are not necessary for the implantation. The application of the products is not depending on the blood type, immunity status etc. of the recipient. The products of CTBA are free of all potentially immunogene elements.

An implantation in dry condition can damage adjacent cells from the bone. This can be prevented if a rehydrogenation in sterile water will be carried out at least for 30 minutes before the implantation.

The allografts can be bathed to add antibiotics in order to eliminate any infection (for example Tobramycin, Vancomycin, etc.).

It is recommended to curette the recipient bone so that the implantation gets in direct contact with the vital (bleeding) bone.

An optimal integration (favoured press-fit) advantages the healing.

Impaction of cancellous allograft, the following points should be noted:

- Close impaction - high primary stability
- Loose impaction - quick incorporation

The application has to be adapted by the performing surgeon to the local requirements.