


CANCELLOUS CYLINDER

Cancellous bone cylinders are preferably used in sports medicine for filling drill channels in cruciate ligament revisions.

Cancellous cylinder
 Origin: Human
 Tissue: Cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Application: Tunnel filling
 Rehydration: Min. 10 minutes



ORDER INFORMATION*

DESCRIPTION	ITEM NUMBER	D [mm]	H [mm]
Cancellous cylinder	ALO423	10	20
	ALO424	10	30
	ALO425	12	20
	ALO426	12	30
	ALO427	14	20
	ALO428	14	30



* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

**Cells + Tissuebank Austria
gemeinnützige GmbH**

Magnesitstraße 1
A-3500 Krems an der Donau

+43 2732-76954-0
+43 2732-76954-60

office@ctba.at
www.ctba.at