

# CANCELLOUS GRANULES IN THE APPLICATOR

JOINTLY  
BUILDING BRIDGES

**C+TBA**  
austria

## Cancelous cubes

Origin: Human  
Tissue: Cancelous bone  
Processing: Allotec® purification procedure  
Inactivation: Min. SAL10<sup>-6</sup> for viruses and bacteria  
Sterilisation: Gamma irradiation  
Application: Bone void filler  
Rehydration: Min. 10 minutes

The applicator of C+TBA is a special form of primary packaging for cancellous bone granules, which simplifies both the rehydration with a physiological medium as well as the application of the granules into the defect zone.



## ORDER INFORMATION\*

DESCRIPTION	ITEM NUMBER	VOLUME [cc]
Cancelous granules in the applicator	ALO360	7
	ALO361	15
	ALO362	30



\* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

**Cells + Tissuebank Austria  
gemeinnützige GmbH**

Magnesitstraße 1  
A-3500 Krems an der Donau

+43 2732-76954-0  
+43 2732-76954-60

office@ctba.at  
www.ctba.at