


# CANCELLOUS RING

The cancellous ring was specially developed to simplify the reconstruction of the glenoid in a total shoulder prosthesis.

## cancellous ring

Origin: Human  
 Tissue: Cancellouse bone  
 Processing: Allotec® purification procedure  
 Inactivation: Min. SAL10<sup>6</sup> for viruses and bacteria  
 Sterilisation: Gamma irradiaton  
 Application: Remodelling of the glenoid  
 in case of shoulder endoprosthesis  
 Rehydration: Min. 10 minutes



## ORDER INFORMATION\*

DESCRIPTION	iD [mm]	ITEM NUMBER	D [mm]	H [mm]
Cancellous ring	1,5	ALO431	26	10
		ALO433	32	10
		ALO432	26	20
		ALO434	32	20
	7,7	ALO436	26	10
		ALO437	32	10
		ALO435	26	20
		ALO430	32	20



\* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.