



# WEDGES



Wedges are preshaped cancellous or corticocancellous bone grafts, mainly used in corrective osteotomy. C+TBA provides a wide range of wedges (cancellous or cortico-cancellous) with different angles and sizes to precisely address the indication and the patient's individual anatomic preconditions.

**Cancellous wedge**  
 Origin: Human  
 Tissue: Cancellous bone  
 Processing: Allotec® purification procedure  
 Inactivation: Min. SAL10<sup>-6</sup> for viruses and bacteria  
 Sterilisation: Gamma irradiation  
 Application: Corrective osteotomy  
 Rehydration: Min. 10 minutes

**Cortico-cancellous wedge**  
 Origin: Human  
 Tissue: Cortical and cancellous bone  
 Processing: Allotec® purification procedure  
 Inactivation: Min. SAL10<sup>-6</sup> for viruses and bacteria  
 Sterilisation: Gamma irradiation  
 Application: Corrective osteotomy  
 Rehydration: Min. 10 minutes



Cancellous wedge



## ORDER INFORMATION\*

DESCRIPTION	A	ITEM NUMBER	S	D [mm]	H [mm]
Cancellous wedge	7°	ALO462	S	<45	5,0
		ALO460	L	≥45	7,0
	10°	ALO465	S	<45	7,0
		ALO463	L	≥45	10,0
	13°	ALO468	S	<45	10,0
		ALO466	L	≥45	13,0
	16°	ALO470	S	<45	13,0
		ALO469	L	≥45	16,0
Cortico-cancellous wedge	15°	ALO410	-	n.a.	10,0



Cortico-cancellous wedge

\* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

**Cells + Tissuebank Austria  
 gemeinnützige GmbH**

Magnesitstraße 1  
 A-3500 Krems an der Donau

+43 2732-76954-0  
 +43 2732-76954-60

office@ctba.at  
 www.ctba.at